


Specifications

1	Paediatric - Airway Balloon Catheter	It should have a non-compliant Balloon with puncture resistance, the design maintains consistent diameter through multiple dilations. It should have a working length of 42.5 CMS. It should have a flexible shaft with stylet optimised for airway anatomy. Airway Balloons should have atraumatic balloon & stylet tip to crosses narrow stenosis. Airway Balloon should have stylet to help enhanced support & manoeuvrability. Airway Balloon should have malleable shaft can be custom - shaped for easy access & direct visualization. Airway Balloon should reach its stated size at minimum of 12 atm to 16 atm and should withstand 17 atm pressures as per sizes of the Airway balloons. Balloon Catheter should be approved by credible authorities like DGCI, for further safety of patients. Size- 5mm x 24 mm, 7mm x 24mm, 9mm x 24mm	Size	Quantity	30000/- Per Piece
			5mm x 24mm	1 piece	
			7mm x 24mm	2 pieces	
			9mm x 24mm	1 piece	
2	Adult - Airway Balloon Catheter	It should have a non-compliant Balloon with puncture resistance, the design maintains consistent diameter through multiple dilations. It should have a working length of 42.5 CMS. It should have a flexible shaft with stylet optimised for airway anatomy. Airway Balloons should have atraumatic balloon & stylet tip to crosses narrow stenosis. Airway Balloon should have stylet to help enhanced support & manoeuvrability. Airway Balloon should have malleable shaft can be custom - shaped for easy access & direct visualization. Airway Balloon should reach its stated size at minimum 12 atm to 16 atm and should withstand 17 atm pressures as per sizes of the Airway balloons. Balloon Catheter should be approved by credible authorities like DGCI, for further safety of patients.	Size	Quantity	30000/- per piece
			10mm x 40mm	2 Piece	
			12mm x 40mm	2 Piece	
			14mm x 40mm	1 Piece	
			16mm x 40mm	1 Piece	
3	Balloon inflation device	It should have an ergonomic piston handle for safe and controlled inflation/deflation of Airway balloon catheter. It should have clear grading of different pressure levels. It should have high pressure connector tubing to connect to the balloon catheter. It should be compatible with different Airway balloon catheters. It should have clear high pressure cylinder barrel to indicate the level of saline/ fluid inside the inflation device. It should have a piston locking system to prevent accidental inflation/deflation of the Airway balloon catheter. It should be approved by credible authorities for further safety of patients.	01		5000/-


 27/12/2026
 Dr. Anshika Agrawal
 Associate Professor
 Department of Otorhinolaryngology
 AIIMS, Gorakhpur